

## Balmed - Bathtub for balneological baths

### Product short description:



**Balmed** bathtub is a modern bathtub specially designed and intended for balneological treatments - carbonic acid baths or brine baths as well as pearl therapeutic baths (Air 64 option) and bathing using light effects (Chromo option).

### Product gallery:



### Product description:

## Balmed - Bathtub for balneological baths

Bathtub for balneological baths BALMED is specially designed for balneological treatments – carbonic acid or brine baths, also pearl baths (AIR option) or baths with light effects (CHROMO option).

The BALMED bathtub has a quick filling system from the bottom of the basin, which minimizes the release of CO<sub>2</sub> from saturation water during the preparation of a carbonic acid bath.

### Treatment options

Balneotherapy: chloride-sodium baths (brine), sulfide-hydrogen, thermal

Carbonic acid baths using a saturator

Pearl baths - available with AIR64 option

**Features:**

Basin available in two colours: white or calypso green

Basin made of acrylic glass and polyester resins

Made of materials resistant to the influence of chemical and physical factors, which allows for reliable and long-term use in conditions necessary for balneological baths using, e.g. brines, sulfur water

Large diameter of drain for fast emptying of the bathtub

Water installation (except for the connections, valves, thermomixer) is made of PVC which ensures its reliability

Treatment time setting - only with AIR option

CO2-saturated water - delivered via filling nozzles

Water, brine or sulphide water - delivered via spout

Up to 3 types of media: water saturated with CO2, warm or cold brine, warm or cold sulfide water

Control panel (only with option AIR)

Shower head

**Additional options**

**AIR 64**

64 aeration nozzles to pearl baths, air heating system, independent air pump (not possible with mud baths)

**CHROMO 12**

12 single light sources

**CHROMO 24**

24 single light sources

**Footrest**

support for feet

**Available versions:**

| Model  | Number of valves | Warm water | Cold water | Shower head | Medium 1 | Medium 2 | Medium 3 |
|--------|------------------|------------|------------|-------------|----------|----------|----------|
| Balmed |                  |            |            |             |          |          |          |

|                       |   |  |  |  |  |  |  |
|-----------------------|---|--|--|--|--|--|--|
| ECO / STANDARD        |   |  |  |  |  |  |  |
| Balmed ECO / STANDARD | 4 |  |  |  |  |  |  |
| Balmed ECO / STANDARD | 5 |  |  |  |  |  |  |
| Balmed ECO / STANDARD | 6 |  |  |  |  |  |  |

### Technical data:

|   | ECO       | STANDARD  |
|---|-----------|-----------|
| Min. capacity [l]:  | 140       | 150       |
| Capacity (up to overflow) [l]:                            | 240       | 300       |
| Filling* (up to overflow)/emptying** (from overflow) [s]: | 240 / 135 | 240 / 135 |
| Number of filling nozzles:                                | 2         | 2         |
| Power supply (CHROMO/AIR) [V/Hz]:                         | 230/50    | 230/50    |
| Current consumption [A]:                                  | 5         | 5         |

Before making a purchase, read the installation sheet (see the Attachments tab) and be sure of the transport route of the device, taking into account its dimensions and weight.

### Product codes:

Reference: -  
 EAN13: -  
 UPC: -

### Product features:

### Product attributes:

Hydrotherapy - Color: White, Calypso  
 Bathtubs - Chromo: No, 12 single light sources, 24 single light sources  
 Bathtubs - Balmed Models: Balmed ECO, Balmed STANDARD  
 Bathtubs - Balmed Valves: 3, 4, 5, 6  
 Bathtubs - AIR 64: No, Yes



Meden-Inmed Sp. z o.o, Wenedów 2, 75-847 Koszalin | [shop@meden.com.pl](mailto:shop@meden.com.pl) | phone: +48 94 344 90 54

---