

Aquaray Balneo - Bathtub for balneological baths



Product short description:

The **Aquaray Balneo** bathtub is a modern therapeutic bathtub characterized by an advanced structure, designed and intended for hydrotherapy treatments.

Product gallery:



Product description:

Aquaray Balneo - Bathtub for balneotherapy baths

The **Aquaray** bathtub is a modern, therapeutic bathtub with an advanced structure, designed for hydrotherapy treatments. Its main function includes an automatic zoned hydromassage. The hydromassage jets in the **Aquaray** bathtub operate in five, three, or a single zone, depending on the design.

Features:

- Basin available in two colours: white or calypso green
- Material: gelcoat (fiberglass mat reinforced)
- Water usage: minimum 160 L
- Hydromassage: 24 jets operating simultaneously, with adjustable intensity
- Pearl massage: 98 aeration apertures (14x7 configuration)
- Automatic functions: Disinfection, filling, emptying, and rinsing of hydromassage ducts
- Quick-fill system: Located at the bottom of the basin
- User interface: Easily accessible color touch panel
- Water temperature indicator: Displays the temperature of the pouring water during filling; otherwise, shows the temperature of the water in the bathtub
- Sensors: Two water level sensors
- Ergonomic design: Rails for palm support, specially shaped sections for arms and legs to ensure proper body positioning during treatments
- Additional features: Foot support and showerhead

Treatment options:

- Bath with hydromassage
- Brine bath
- Carbonic acid bath (using a saturator)
- Beer bath
- Pearl bath (available with **AIR98** option)

Additional options:

- AIR 98: 98 aeration apertures for pearl massage (14x7 apertures)
- CHROMO 12: 12 single light sources
- CHROMO 24: 24 single light sources
- CHROMO 36: 36 single light sources
- AUDIO: Sound input from an external device such as Hi-Fi, MP3, etc.
- OZONE: ozone generator
- THERMOSTAT: thermostatic battery (replaces the manual hot and cold water filling valves)

Available versions:

Number of valves	Warm water	Cold water	Shower head	CO2	Warm saline	Cold saline
3						
4						
5						
6						

Technical data:

Min. bath level capacity [l]:	160
Capacity up to overflow [l]:	320
Filling (up to min. bath level)/emptying [min]:	4 / 2,5

Current consumption [a]:	13
Power supply [V/Hz]:	230/50
Dimensions (L x W x H) [cm]:	225 x 88 x 91
Weight [kg]:	160
Max. patient weight [kg]::	135
Protection class, applied part:	Class I, Type B
Zone hydromassage pump - pressure range [bar]:	1,5 - 2,2

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Hydrotherapy - Color: White, Calypso
 Bathtubs - Chromo: No, 12 single light sources,
 24 single light sources, 36 single light sources
 Bathtubs - Audio: No, Yes
 Bathtubs - AIR 98: No, Yes
 Bathtubs - Ozone: No, Yes
 Bathtubs - Thermostat: No, Yes