

Fiji - Pacific examination and treatment tables



Product short description:

Fiji - examination and treatment table is a modern **5 and 7-section table** that features **electric height adjustment, a heated table surface and LED lighting - (optional)**. It is reliable, aesthetically pleasing, and comfortable to use. The table's modern design will enhance the prestige of your practice and make your clients feel valued.

Product gallery:



Product description:

Fiji - Pacific examination and treatment table

This modern 5 and 7-section examination and treatment table features an ergonomic design that ensures therapist comfort. The wide lying surface ensures patient safety and comfort.

Features:

- 5- or 7-section examination and treatment table with electric height adjustment by the means of a remote

(controller).

- Head section adjustment by the means of gas spring.
- Head rest with shaped face hole with crescent plug.
- Foldable armrests lowered till 19 cm (only in 7 section table).
- Adjustment of: back, thighs, shanks sections by means of electric actuators.
- Multi-layered and very soft padding of 15 cm thickness with upholstery in many colours.
- Trendelenburg/reverse Trendelenburg position (15°).
- 2 wheels.

Additional options:

- An easy-to-use heating panel with two temperature settings for the padding of the **Fiji examination and treatment table**.
- **LED lighting**, which allows you to choose your preferred lighting color to match the interior of your medical office.
- **The central mechanism lifts the transport wheels.**

The remote control is designed to be simple and easy to use

Fiji - examination and treatment table is equipped with a 10-button remote control that allows quick repositioning of the chair:

- Backrest adjustment
- Thigh section adjustment
- Shank section adjustment
- Height adjustment of the table
- Tilt table (chair function)

Fiji is a luxury examination and treatment table covered with unique upholstery that is:

- Stretch-resistant
- Extremely soft and velvety
- Easy to clean
- Covered with **Silverguard®** for fungal and bacterial resistance.
- Resists many liquids, including **alcohol, oil, blood, and urine.**
- Low temperature resistant
- UV resistant
- Abrasion resistant

Technical data:

Numebr of sections	5 or 7
Height adjustment [cm]:	62-10259-99
Height adjustment of armrests [cm]:	0-19
Headrest section angle adjustment [°]:	0/-500/+38
Chest section angle adjustment [°]:	0/+42

Thigh section angle adjustment [°]:	0/+25
Shank section angle adjustment [°]:	0/+18
Trendelenburg / anti-Trendelenburg position [°]:	0/+150/-15
Upholstery thickness [cm]:	9 or 15
Dimensions (L x W) [cm]:	212 x 80
Weight [kg]:	150
Max. lifting capacity [kg]:	150
Power supply [V/Hz]:	220-240 / 50-60
Medical device:	

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Type of table: Examination table

Type of table: Treatment table

Number of sections: 5

Number of sections: 7

Wheeling system: F0

Length [cm]: 212

Width [cm]: 80

Weight [kg]: 150

Max. lifting capacity [kg]: 150

Power supply [V / Hz]: 230 / 60

Warranty: 24 months

Medical device: Yes

Product attributes:

Tables - Fiji / Hawaii: 5 - section treatment table, 7 - section treatment table

Tables - Additional options: -, Heated table surface + LED Lightning