

## Elemental Ice E-S2 - Massage table

### Product short description:

**The Elemental Ice E-S2** is a two-section massage and rehabilitation table with electric height adjustment via remote control. It comes equipped with a basic wheel system consisting of 2 small wheels and 2 non-slip rubber feet with a height adjustment of up to 1 cm, allowing the table to be precisely positioned.



### Product gallery:





## Product description:

### Elemental Ice E-S2 - Massage table

The **Elemental Ice E-S2** massage table is designed for performing manual massage treatments to treat and alleviate traumas or disabilities by restoring or improving the function of massaged tissues, organs and systems.

The clinical benefits of massage tables are based on the fact that it is possible for the patient to take the right position for the massage treatment, the possibility to adjust the height and/or individual sections of the table, which facilitates the work of the rehabilitator/physiotherapist and increases the comfort of the patient during the massage treatment.

**Elemental Ice E-S2** - Massage table by **Meden-Inmed** is appreciated by therapists all over the world.

## Features:

2-section massage table.

Adjustable head, leg section with gas springs.

Round face hole in the head section with 2 removable face hole covers (crescent, round).

2 side brackets for the fixation of stabilization belts

Electric height adjustment: 48-96 cm via hand-operated remote controller (with authorization) or foot controller to choose in standard

High/low height adjustment foot bar accessible from each side of the table - option (replaces remote controller).

Maximum lifting capacity up to 180kg.

Maximum working load up to 250kg.

Locking safety system on the remote controller.

## Types of base

Massage tables from Elemental Series are equipped with two types of wheeling systems. You can choose between 2 small wheels with rubber feet or add the integrated F2 mechanism, which adds two additional wheels for increased maneuverability and seamless transport.

### F0 - Wheeling system

Constructed with durable materials, the Elemental Table is equipped with 2 legs featuring rubber feet and 2 small wheels, facilitating effortless relocation by lifting from one side. This feature-rich system is ideally suited for users who only occasionally need to move their Elemental Table, ensuring both stability and convenience in its use.

### F0+F2 - Wheeling system

The innovative design of the massage table ensures stability and ease of movement thanks to the combination of two small wheels and two sturdy feet. The integrated F2 mechanism, consisting of two casters, enhances maneuverability and allows seamless transport and positioning of the table in various therapy settings.

Elemental Series

## A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Elemental Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

The Elemental Series massage table features:

Two-layer padding (3 cm foam 90 kg/cm<sup>3</sup> + 1 cm foam 43 kg/cm<sup>3</sup>).

Stitched upholstery

A frame combining wood and steel.

High-quality, soft upholstery in a wide variety of colors

Upholstery in Elemental Series

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

## Temperature resistant

The two-layer upholstery of the **Elemental Ice E-S2** massage table is resistant to high temperatures.

## Liquid resistant

The top material of the **Elemental Ice E-S2** massage table is resistant to liquids (water, blood, brine).

## Comfortable surface

The **Elemental Ice E-S2** massage table is covered with a soft, flexible and comfortable coating.

## Scratch resistant

The two-layer upholstery of the **Elemental Ice E-S2** features a scratch-resistant structure.



## Product codes:

Reference: -

EAN13: -

UPC: -

## Product features:

Type of table: Rehabilitation table

Number of sections: 2

Wheeling system: F0

Wheeling system: F2

Pivot: No

Length [cm]: 202

Width [cm]: 69

Height [cm]: 48-96

Weight [kg]: 81

Max. lifting capacity [kg]: 180

Max. working load [kg]: 250

Power supply [V / Hz]: 230 / 60

Warranty: 24 months

Medical device: Yes

## Product attributes:

Tables - Color: OM006 - White, OM240 - Beige,

OM315 - Blue, OM221 - Dark blue, OM314 -

Grey, OM076 - Graphite, OM173 - Brown,

OM225 - Black, OM469 - Light green, OM157 -

Green, OM053 - Red, OM217 - Burgundy,

OM095 - Mostaza, OM075 - Orange

Tables - Upholstery thickness: 40 mm, 50 mm

Tables - Additional controls: No, Foot Remote

Control, Frame