



## Vertimo Classic - Tilting table

### Product short description:

Vertimo - Tilting table is a series of tables designed to help patients stand upright. These tables come with electric height and angle adjustment features. The tables feature a 2-layer upholstery made of approved materials, available in many color options.

### Product gallery:

### Product description:

## Vertimo Classic - Tilting table

**Vertimo Classic** is a stable tilting table with electric angle adjustment. Its one-section design and fixed height of 53 cm facilitate patient transfer from a wheelchair or bed.

### Features:

- 1-section table.
- Electric actuator for tilting adjustment.
- Manual remote control for the adjustment of tilting angle with access authorization - with a special key placed in the remote control.
- Brackets for the attachment of fastening belts (standard accessories).
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable steel frame.
- 4 large separately lockable castors (rolling and rotation lock).
- Emergency power supply battery.

### Technical data:

Dimensions (Length x Width) [cm]:	190 x 69
Verticalization angle adjustment [°]:	0-87 (+3)
Height [cm]:	53

Max. working load [kg]:	200
Weight [kg]:	75
Medical device:	

**Product codes:**

Reference: 01-STO-2020-001#

EAN13: -

UPC: -

**Product features:**

Type of table: Tilting table

Number of sections: 1

Wheeling system: F4

Pivot: No

Length [cm]: 190

Width [cm]: 69

Height [cm]: 53

Weight [kg]: 75

Max. working load [kg]: 200

Warranty: 24 months

Medical device: Yes

**Product attributes:**