

Vertimo Hi-LO - Tilting table



Product short description:

Vertimo is a series of patient tilting tables with electric height and angle adjustment. The tables are equipped with a 2-layer upholstery made of approved materials and are available in a wide range of colour options.

Product gallery:

Product description:

Vertimo Hi-Lo - Tilting table

The **Vertimo Hi-Lo** is a very popular patient tilting table in the field, offering the following features **easy pivot angle** and **height adjustment** for working with neurological patients.

Features:

- 1-section table.,
- 2 electric actuators – one ensures smooth tilting; the other one adjusts height.
- Electric height and tilting adjustment with manual remote control with the function of access authorization - authorization is done by a special key placed in the remote control.
- Two-layer padding with upholstery in many colours.
- Powder coated, very stable steel frame.
- 2 separately adjustable footrests.
- Brackets for the attachment of fastening belts (standard accessories).
- 4 large separately lockable castors (rolling and rotation lock).
- Tilt angle indicator.
- Emergency power supply battery (allows table lowering to min. height and return to horizontal position at power failures).

Technical data:

Dimensions (Length x Width) [cm]:	190 x 69
-----------------------------------	----------

Verticalization angle adjustment [°]:	0-87 (+3)
Height adjustment [cm]:	53-103
Max. working load [kg]:	200 (150 in vertical angle)
Weight [kg]:	115
Medical device:	

Product codes:

Reference: 01-STO-2020-001#DOL18-R9002DS-H40NS-MED-EU

EAN13: -

UPC: -

Product features:

Type of table: Tilting table

Number of sections: 1

Wheeling system: F4

Pivot: No

Length [cm]: 190

Width [cm]: 69

Height [cm]: 53-103

Weight [kg]: 115

Max. working load [kg]: 200

Warranty: 24 months

Medical device: Yes

Product attributes: