

Prestige E-S2 - Massage table



Product short description:

Prestige E-S2 massage table is equipped with the innovative **Hallotronic®** system and **Personal Authorization System** for comfortable and convenient work with the patient.

The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with non-slip rubber feet that can be adjusted by up to 1 cm.

Product gallery:

Product description:

Prestige E-S2 - Massage table

Prestige E-S2 massage table is equipped with the innovative **Hallotronic®** system and **Personal Authorization Safety System (PASS)** for comfortable and convenient work with the patient. The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with rubber feet that can be adjusted by up to 1 cm.

Prestige E-S2 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

Prestige Series

A unique experience for your patients

For the patient, comfort during the treatment is very important. With the Prestige Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

The Prestige E-S2 massage table features:

2-section massage table with a long head section

Adjustable head section from 0° to +85° with 2 gas springs.

Face hole in the head section with 2 removable face hole covers (round, crescent).

Optional 2 side brackets for the fixation of stabilization belts and elastic bands or 2 retractable side railings for fall prevention.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Personal Authorization Safety System with 2 safety magnetic keys.

HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg

Two ways of frame painting: „white” or „twintone”.

Types of base

You can choose between F0 and F4 wheeling system. F0 system composed of 2 legs and 2 small wheels allows for moving of Prestige Table by lifting it from one side. The system is recommended for users who relocate their Prestige Table less frequently. F4 system composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

F0 - Wheeling system

Composed of 2 legs with rubber feet and 2 small wheels for moving of Prestige Table by lifting it from one side. It is recommended for users who relocate their Prestige Table less frequently.

F4 - Wheeling system

Composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

Technologies

Personal Authorization Safety System

Tables have an integrated electronic control unit allowing for height adjustment and adjustment of the postural drainage option with an indication of stand-by mode (green light) and Personal Authorization Safety System (PASS). Each Prestige Table is equipped with two magnetic safety keys (PASS Keys) in order to prevent unauthorized use of a table which may lead to severe accidents and compensation claims.

Hallotronic® Technology

We care about your working comfort! That's why we invented the unique HalloTronic® System that enables high/low height adjustment via hallotrons without the use of mechanical switches. Thanks to HalloTronic® System you can choose the right height of your Prestige Table without effort and noise.

Prestige - Massage tables are widely used in:

Sanatorium

Massage tables from the Prestige Series are a top-quality choice for sanatoriums, offering a comprehensive range of solutions designed to meet the needs of even the most demanding customers. Built with advanced technology and precision craftsmanship, these tables fully **support the recovery and relaxation process** of the

client.

Rehabilitation clinics

Prestige massage tables are unparalleled in their class, providing comprehensive solutions for rehabilitation centres that are critical to the recovery process and **improving patients' quality of life**. These tables offer not only comfort but also **effective therapy**, thanks to their advanced technology and precise workmanship.

Upholstery in Prestige Series

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

Easy to clean

We recommend wiping **Prestige E-S2** with a neutral soap solution using a cloth in smooth circular movements.

Liquid resistant

The top material of the **Prestige E-S2** massage table is resistant to liquids (water, blood, brine).

Comfortable surface

The **Prestige E-S2** massage table is covered with a soft, flexible and comfortable coating.

Scratch resistant

The top coating of the **Prestige E-S2** features a scratch-resistant structure.

Upholstery thickness available in two versions:

40

mm

80

mm

40
80
120
160
200
240
280
320
360
400
440
480
520
560
600
640
680
720
760
800
840
880
920
960
1000
1040
1080
1120
1160
1200
1240
1280
1320
1360
1400
1440
1480
1520
1560
1600
1640
1680
1720
1760
1800
1840
1880
1920
1960
2000
2040
2080
2120
2160
2200
2240
2280
2320
2360
2400
2440
2480
2520
2560
2600
2640
2680
2720
2760
2800
2840
2880
2920
2960
3000
3040
3080
3120
3160
3200
3240
3280
3320
3360
3400
3440
3480
3520
3560
3600
3640
3680
3720
3760
3800
3840
3880
3920
3960
4000
4040
4080
4120
4160
4200
4240
4280
4320
4360
4400
4440
4480
4520
4560
4600
4640
4680
4720
4760
4800
4840
4880
4920
4960
5000
5040
5080
5120
5160
5200
5240
5280
5320
5360
5400
5440
5480
5520
5560
5600
5640
5680
5720
5760
5800
5840
5880
5920
5960
6000
6040
6080
6120
6160
6200
6240
6280
6320
6360
6400
6440
6480
6520
6560
6600
6640
6680
6720
6760
6800
6840
6880
6920
6960
7000
7040
7080
7120
7160
7200
7240
7280
7320
7360
7400
7440
7480
7520
7560
7600
7640
7680
7720
7760
7800
7840
7880
7920
7960
8000
8040
8080
8120
8160
8200
8240
8280
8320
8360
8400
8440
8480
8520
8560
8600
8640
8680
8720
8760
8800
8840
8880
8920
8960
9000
9040
9080
9120
9160
9200
9240
9280
9320
9360
9400
9440
9480
9520
9560
9600
9640
9680
9720
9760
9800
9840
9880
9920
9960
10000

Additional accessories

Wedges

Bolsters

Paper holder

Brackets/Railings

Stabilization Belt

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Type of table: Massage table

Number of sections: 2

Wheeling system: F0

Wheeling system: F4

Pivot: No

Length [cm]: 200

Width [cm]: 69

Height [cm]: 47-95

Weight [kg]: 115

Max. working load [kg]: 300

Warranty: 24 months

Medical device: Yes

Product attributes: