

## Prestige M-P3 - Massage table



### Product short description:

**Prestige M-P3** massage table is equipped with the innovative **Hallotronic®** system and **Personal Authorization Safety System (PASS)** for comfortable and convenient work with the patient.

The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with non-slip rubber feet that can be adjusted by up to 1 cm.

### Product gallery:

### Product description:

## Prestige M-P3 - Massage table

**Prestige M-P3** massage and rehabilitation table is equipped with the innovative **Hallotronic®** system and **Personal Authorization Safety System (PASS)** for comfortable and convenient work with the patient. The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with rubber feet that can be adjusted by up to 1 cm.

Prestige M-P3 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

### Prestige Series

For the patient, comfort during the treatment is very important. With the Prestige Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

### The Prestige M-P3 massage table features:

3-section massage and rehabilitation table

Pivot option - postural drainage with an electric actuator operated by 2 switches from both sides of the table.

Adjustable head section from -70° to +40° with a gas spring.

Face hole in the head section with 2 removable face hole covers (round, crescent).

Adjustable leg section from -10° to 80° with 2 gas springs.

Release handle of the gas spring system in the leg section accessible from both sides of the table.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Personal Authorization Safety System with 2 safety magnetic keys.

HalloTronic® System for electric high/low height adjustment from 51 to 99 cm with a lifting capacity up to 225 kg

Traction support sleeves.

High/low height adjustment foot bar accessible from each side of the table.

Two ways of frame painting: „white” or „twintone”.

## **Types of base**

You can choose between F0 and F4 wheeling system. F0 system composed of 2 legs and 2 small wheels allows for moving of Prestige Table by lifting it from one side. The system is recommended for users who relocate their Prestige Table less frequently. F4 system composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

### **F0 - Wheeling system**

Composed of 2 legs with rubber feet and 2 small wheels for moving of Prestige Table by lifting it from one side. It is recommended for users who relocate their Prestige Table less frequently.

### **F4 - Wheeling system**

Composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

## **Technologies**

### **Personal Authorization Safety System**

Tables have an integrated electronic control unit allowing for height adjustment and adjustment of the postural drainage option with an indication of stand-by mode (green light) and Personal Authorization Safety System (PASS). Each Prestige Table is equipped with two magnetic safety keys (PASS Keys) in order to prevent unauthorized use of a table which may lead to severe accidents and compensation claims.

### **Hallotronic® Technology**

We care about your working comfort! That's why we invented the unique HalloTronic® System that enables high/low height adjustment via hallotrons without the use of mechanical switches. Thanks to HalloTronic® System you can choose the right height of your Prestige Table without effort and noise.

### **Prestige - Massage tables are widely used in:**

## Sanatorium

Massage tables from the Prestige Series are a top-quality choice for sanatoriums, offering a comprehensive range of solutions designed to meet the needs of even the most demanding customers. Built with advanced technology and precision craftsmanship, these tables fully **support the recovery and relaxation process** of the client.

## Rehabilitation clinics

Prestige massage tables are unparalleled in their class, providing comprehensive solutions for rehabilitation centres that are critical to the recovery process and **improving patients' quality of life**. These tables offer not only comfort but also **effective therapy**, thanks to their advanced technology and precise workmanship.

## Upholstery in Prestige Series

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

## Easy to clean

We recommend wiping **Prestige M-P3** with a neutral soap solution using a cloth in smooth circular movements.

## Liquid resistant

The top material of the **Prestige M-P3** massage table is resistant to liquids (water, blood, brine).

## Comfortable surface

The **Prestige M-P3** massage table is covered with a soft, flexible and comfortable coating.

## Scratch resistant

The two-layer upholstery of the **Prestige M-P3** features a scratch-resistant structure.

## Upholstery thickness available in two versions:

**40**

mm

**80**

mm

- Standard Upholstery
- Premium Upholstery
- Omega Upholstery

## Technical specification:

**Upholstery: eco leather** Highest quality and functionality

Liquid (water, blood) resistant High temp. resistant Soft, elastic and comfortable Scratch resistant Easy to clean

**00 - Ecu**

**01 - Beige**

**02 - Blue**

**03 - Navy Blue**

**04 - Terracotta**

**05 - Brown**

**06 - Grey**

**07 - Graphite**

**08 - Black**

**09 - Red**

**11 - Honey**

**12 - Sea Green**

**13 - Orange**

**14 - Lime**

**15 - Light Blue**

### **Technical specification:**

**Coating: 100% PolyurethanePhthalate - free**

6 times more durable Low temp. resistant UV Resistant Chemical resistant Abrasion resistant

**White - Dolce 6701**

**Ecu - Dolce 6703**

**Mango - Dolce 6731**

**Grass - Dolce 6722**

**Jet - Dolce 6724**

**Quetzal - Dolce 6734**

**Tangelo - Dolce 6718**

**Wine - Dolce 6738**

**Artesian - Dolce 6726**

**Graphite - Dolce 6739**

**Black - Dolce 6728**

**Chocolate - Dolce 6713**

### **Technical specification:**

**Surface: 100% PolyurethaneBacking: 100% Polyester**

Resistance to hydrolis (8 years) No chemical solvents High scratch resistance Respects the environment  
Durability and performance

**OM006 - White**

**OM240 - Beige**

**OM315 - Sky**

**OM221 - Blue**

**OM314 - Cement**

**OM076 - Olive Grey**

**OM173 - Dark Brown**

**OM225 - Black**

**OM469 - Lima**

**OM157 - Forest**

**OM053 - Red**

**OM217 - Bordeaux**

### **Additional accessories**

**Wedges**

**Bolsters**

**Paper holder**

**Brackets/Railings**

**Stabilization Belt**

### **Product codes:**

Reference: -

EAN13: -

UPC: -

### **Product features:**

Type of table: Massage table

Number of sections: 3

Wheeling system: F4

Pivot: Yes

Length [cm]: 205

Width [cm]: 69

Height [cm]: 51-99

Weight [kg]: 110

Max. lifting capacity [kg]: 225

Max. working load [kg]: 300

### **Product attributes:**

Tables - Color: 00 - Ecru - Standard, 01 - Beige - Standard, 02 - Blue- Standard, 03 - Navy Blue - Standard, 04 - Teracota - Standard, 05 - Brown - Standard, 06 - Grey - Standard, 07 - Graphite - Standard, 08 - Black - Standard, 09 - Red - Standard, 11 - Honey - Standard, 12 - Marine - Standard, 13 - Orange - Standard, 14 - Lime - Standard, 15 - Navy - Standard, DOL-6701 White - Premium, DOL-6703 Ecru - Premium, DOL-6731 Mango - Premium,

Power supply [V / Hz]: 230 / 60

Warranty: 24 months

Medical device: No

DOL-6722 Grass - Premium, DOL-6724 Jet - Premium, DOL-6734 Quetzal - Premium, DOL-6726 Artesian - Premium, DOL-6739 Graphite - Premium , DOL-6728 Black - Premium, DOL-6718 Tangelo - Premium, DOL-6738 Wine - Premium, DOL-6713 Chocolate - Premium

Tables - Upholstery thickness: 40 mm, 80 mm