

MoHiCan II - Hygiene and toilet chair

Product short description:

MoHiCan II hygiene and toilet chair is designed to assist people with limited mobility to move and perform daily hygiene and toilet activities.

[View product presentation](#)



Product gallery:





Product description:

MoHiCan II - Hygiene and toilet chair

MoHiCan II mobile hygienic chair with adjustable position is used as an aid for transportation and performance of hygienic and toileting procedures in institutional and home care of elderly people, people with paresis of lower limbs or people with other disabilities that make it difficult to perform hygienic and toileting procedures independently.

This mobile hygienic chair has electric height and patient positioning adjustments and offers protection against water, dirt and hair flow after a patient's haircut or shave. The positions of the various parts of the chair are adjustable via a handheld remote control, ensuring optimal patient positioning.

Mohican II facilitates access to the shower and assists with transportation, hygiene, and toileting for the elderly, persons with paralysis of the lower extremities, or other disabilities that make independent hygiene and toileting difficult.

Features

Electric height adjustment and forward and backward tilt of the seat,

Reclining armrests for comfort and easier access to the patient,
Adjustable headrest that provides head stability for greater patient comfort,
4 double casters with brake for safety and convenience of use,
Removable seat makes it easy to perform sanitary tasks,
Footrest for greater patient comfort,
Maximum safe working load up to 180kg,
Lifting range: 47,6 - 100,5,

INTERACTIVE 360° PRODUCT EXPERIENCE

Drag to rotate **Scroll** to zoom 2x **Double click** for full screen

MoHiCan II will be perfect for:

Hospital

MoHiCan II is a valuable addition to any hospital, providing exceptional comfort and hygiene features for patients in need of special care.

Care home

Its hygiene features, including removable seat cushions and washable upholstery, ensure optimal cleanliness in a care home setting.

Private usage

MoHiCan II hygienic chair is also suitable for private usage. The chair's flexibility allows for easy maneuverability, making it convenient for individuals to transfer in and out of the chair.

Sanatorium

MoHiCan II hygienic chair is well-suited for use in sanatoriums. The durability of the chair ensures its suitability for long-term use in high-demand settings such as sanatoriums.

Additional accessories:

Removable container for hygienic and toilet treatment.

Trunk stabilization belt helps keep patient in proper position

Stabilizing leg strap provides additional protection and stabilization

Technical data:

Length [cm]:	89
Width [cm]:	65,8
Handle height [cm]:	102,8 - 126
Lifting height [cm]:	49,5 - 101
Tilt adjustment [degrees]:	-5 / +27
Footrest height from the floor [cm]:	11,2
Turning diameter of the device [cm]:	120
Caster diameter [cm]:	12,5
Weight [kg]:	60
Battery capacity [Ah]:	2,9
Charging time:	approx. 5 hours
Maximum load [kg]:	180

Video presentation:

Product codes:

Reference: 01-KPR-0101-001

EAN13: -

UPC: -

Product features:

Product attributes: