



MoHiCan II Plus - commode and shower chair

Product short description:

MoHiCan II Plus is a mobile, adjustable-position hygiene chair with a load capacity of up to 250 kg. It is designed to help individuals with reduced mobility perform daily personal hygiene and toileting activities.

Product gallery:





Product description:

MoHiCan II Plus - Hygiene and toilet chair

MoHiCan II Plus mobile hygienic chair with adjustable position is used as an aid for transportation and performance of hygienic and toileting procedures in institutional and home care of elderly people, people with paresis of lower limbs or people with other disabilities that make it difficult to perform hygienic and toileting procedures independently.

This mobile hygienic chair has electric height and patient positioning adjustments and offers unique protection against water, dirt and hair flow after a patient's haircut or shave. The positions of the various parts of the chair are adjustable via a handheld remote control, ensuring optimal patient positioning.

Mohican II Plus facilitates access to the shower and assists with transportation, hygiene, and toileting for the elderly, persons with paralysis of the lower extremities, or other disabilities that make independent hygiene and toileting difficult.

Indications for using the MoHiCan II Plus hygiene chair

The MoHiCan II Plus hygiene chair does not have a direct clinical effect. It is useful in aiding staff (assistants) in caring for patients with conditions such as:

MoHiCan II Plus will per
Paralysis,
Disability,
Musculoskeletal changes,
Chronic immobilization,
Parkinson's disease,
Craniocerebral injuries,
Mobility disabilities,
Bedsores,
Guono,

MoHiCan II Plus will perfectly fit in:

Hospital

Stroke

MoHiCan II Plus is a valuable addition to any hospital, providing exceptional comfort and hygiene features for patients in need of special care.

Care home

MoHiCan II Plus has hygiene features such as removable seat cushions and washable upholstery, which ensure optimal cleanliness in a care home setting.



Private usage

MoHiCan II Plus hygienic chair is also suitable for private usage. The chair's flexibility allows for easy maneuverability, making it convenient for individuals to transfer in and out of the chair.

Sanatorium

MoHiCan II Plus hygienic chair is well-suited for use in sanatoriums. The durability of the chair ensures its suitability for long-term use in high-demand settings such as sanatoriums.

Additional accessories:

Removable container for hygienic and toilet treatment.

Trunk stabilization belt helps keep patient in proper position

Stabilizing leg strap proivdes additional protection and stabilization

Technical data:

Length [cm]:	89
Width [cm]:	77,8
Handle height [cm]:	102,8 - 126
Lifting height [cm]:	49 - 101
Tilt adjustment [degrees]:	-5 / +27
Footrest height from the floor [cm]:	11,2
Turning diameter of the device [cm]:	120
Caster diameter [cm]:	12,5
Weight [kg]:	70
Battery capacity [Ah]:	2,9
Charging time:	approx. 5 hours
Maximum load [kg]:	250

Video presentation:

Product codes:

Reference: -

EAN13: -

UPC: -



Product features:	Product attributes:
-------------------	---------------------