

## Prestige M-S5 - Massage table

### Product short description:

Prestige M-S5 massage table is equipped with the innovative **Hallotronic®** system and Personal Authorization System for comfortable and convenient work with the patient.

The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with non-slip rubber feet that can be adjusted by up to 1 cm.

Prestige M-S5 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

### Product gallery:



### Product description:

Prestige M-S5 - Massage table

**Prestige M-S5** massage and rehabilitation table is equipped with the innovative **Hallotronic®** system and

**Personal Authorization System** for comfortable and convenient work with the patient. The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with non-slip rubber feet that can be adjusted by up to 1 cm.

Prestige M-S5 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

Prestige Series

### **A unique experience for your patients**

For the patient, comfort during the treatment is as important as the results. With the Prestige Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

### **The Prestige M-S5 massage table features:**

5-section massage and rehabilitation table

Adjustable head section from -70° to +40° with a gas spring.

Face hole in the head section with 2 removable face hole covers (round, crescent).

Fold-down side sections in the middle section with a semi-automatic blocking/release system.

Release handle of the gas spring system in the leg section accessible from both sides of the table.

Traction support sleeves.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Personal Authorization Safety System with 2 safety magnetic keys.

HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg

Innovative self-locking fold-down armrests in head section (one hand operated) - SmartS4Arms (Smart Slide For Arms).

### **Types of base**

You can choose between F0 and F4 wheeling system. F0 system composed of 2 legs and 2 small wheels allows for moving of Prestige Table by lifting it from one side. The system is recommended for users who relocate their Prestige Table less frequently. F4 system composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

#### **F0 - Wheeling system**

Composed of 2 legs with rubber feet and 2 small wheels for moving of Prestige Table by lifting it from one side. It is recommended for users who relocate their Prestige Table less frequently.

#### **F4 - Wheeling system**

Composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

Technologies

## Personal Authorization Safety System

Tables have an integrated electronic control unit allowing for height adjustment and adjustment of the postural drainage option with an indication of stand-by mode (green light) and Personal Authorization Safety System (PASS). Each Prestige Table is equipped with two magnetic safety keys (PASS Keys) in order to prevent unauthorized use of a table which may lead to severe accidents and compensation claims.

## Hallotronic® Technology

We care about your working comfort! That's why we invented the unique HalloTronic® System that enables high/low height adjustment via hallotrons without the use of mechanical switches. Thanks to HalloTronic® System you can choose the right height of your Prestige Table without effort and noise.

## SmartS4Arms Technology

Innovative system for smooth, screwless adjustment of the headrest sides, based on a sliding and locking mechanism. The automatic position lock prevents any unwanted movement of the headrest sides during therapy.

## Prestige - Massage tables are widely used in:

### Sanatorium

Massage tables from the Prestige Series are a top-quality choice for sanatoriums, offering a comprehensive range of solutions designed to meet the needs of even the most demanding customers. Built with advanced technology and precision craftsmanship, these tables fully **support the recovery and relaxation process** of the client.

### Rehabilitation clinics

Prestige massage tables are unparalleled in their class, providing comprehensive solutions for rehabilitation centres that are critical to the recovery process and **improving patients' quality of life**. These tables offer not only comfort but also **effective therapy**, thanks to their advanced technology and precise workmanship.

## Upholstery in Prestige Series

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

### Temperature resistant

The two-layer upholstery of the **Prestige M-S5** massage table is resistant to high temperatures.

### Liquid resistant

The top material of the **Prestige M-S5** massage table is resistant to liquids (water, blood, brine).

### Comfortable surface

The **Prestige M-S5** massage table is covered with a soft, flexible and comfortable coating.

## Scratch resistant

The two-layer upholstery of the **Prestige M-S5** features a scratch-resistant structure.

## Upholstery thickness available in two versions:

**40**

mm

**80**

mm



## Accessories

### Wedges

### Bolsters

### Paper holder

### Protection mat

### Stabilization Belt

## Product codes:

Reference: -

EAN13: -

UPC: -

## Product features:

Type of table: Massage table

Number of sections: 5

Wheeling system: F4

Pivot: No

Length [cm]: 205

Width [cm]: 69

Height [cm]: 47-95

Weight [kg]: 107

Max. lifting capacity [kg]: 225

Max. working load [kg]: 300

Power supply [V / Hz]: 230 / 60

## Product attributes:

Warranty: 36 months

Medical device: Yes