

Prestige M-S2 - Massage table



Product short description:

Prestige M-S2 massage table is equipped with the innovative **Hallotronic®** system and Personal Authorization System for comfortable and convenient work with the patient.

The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with non-slip rubber feet that can be adjusted by up to 1 cm.

Prestige M-S2 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

Product gallery:

Product description:

Prestige M-S2 - Massage table

Prestige M-S2 massage and rehabilitation table is equipped with the innovative **Hallotronic®** system and **Personal Authorization System** for comfortable and convenient work with the patient. The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with non-slip rubber feet that can be adjusted by up to 1 cm.

Prestige M-S2 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

Prestige Series

A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Prestige Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

The Prestige M-S2 massage table features:

2-section massage and rehabilitation table

Adjustable head section from -70° to +40° with a gas spring.

Face hole in the head section with 2 removable face hole covers (round, crescent).

Adjustable leg section from -10° to 80° with 2 gas springs.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Personal Authorization Safety System with 2 safety magnetic keys.

HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg

Types of base

You can choose between F0 and F4 wheeling system. F0 system composed of 2 legs and 2 small wheels allows for moving of Prestige Table by lifting it from one side. The system is recommended for users who relocate their Prestige Table less frequently. F4 system composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

F0 - Wheeling system

Composed of 2 legs with rubber feet and 2 small wheels for moving of Prestige Table by lifting it from one side. It is recommended for users who relocate their Prestige Table less frequently.

F4 - Wheeling system

Composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

Technologies

Personal Authorization Safety System

Tables have an integrated electronic control unit allowing for height adjustment and adjustment of the postural drainage option with an indication of stand-by mode (green light) and Personal Authorization Safety System (PASS). Each Prestige Table is equipped with two magnetic safety keys (PASS Keys) in order to prevent unauthorized use of a table which may lead to severe accidents and compensation claims.

Hallotronic® Technology

We care about your working comfort! That's why we invented the unique HalloTronic® System that enables high/low height adjustment via hallotrons without the use of mechanical switches. Thanks to HalloTronic® System you can choose the right height of your Prestige Table without effort and noise.

Massage tables available used in:

Massage tables from the Prestige Series are a top-quality choice for sanatoriums, offering a comprehensive range of solutions designed to meet the needs of even the most demanding customers. Built with advanced technology and precision craftsmanship, these tables fully **support the recovery and relaxation process** of the client.

Prestige massage tables are unparalleled in their class, providing comprehensive solutions for rehabilitation centres that are critical to the recovery process and **improving patients' quality of life**. These tables offer not only comfort but also **effective therapy**, thanks to their advanced technology and precise workmanship.

Prestige - Massage tables are widely used in:

Sanatorium

Massage tables from the Prestige Series are a top-quality choice for sanatoriums, offering a comprehensive range of solutions designed to meet the needs of even the most demanding customers. Built with advanced technology and precision craftsmanship, these tables fully **support the recovery and relaxation process** of the client.

Rehabilitation clinics

Prestige massage tables are unparalleled in their class, providing comprehensive solutions for rehabilitation centres that are critical to the recovery process and **improving patients' quality of life**. These tables offer not only comfort but also **effective therapy**, thanks to their advanced technology and precise workmanship.

Upholstery in Prestige Series

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

Temperature resistant

The two-layer upholstery of the **Prestige E-S2** massage table is resistant to high temperatures.

Liquid resistant

The top material of the **Prestige E-S2** massage table is resistant to liquids (water, blood, brine).

Comfortable surface

The **Prestige E-S2** massage table is covered with a soft, flexible and comfortable coating.

Scratch resistant

The two-layer upholstery of the **Prestige E-S2** features a scratch-resistant structure.

Upholstery thickness available in two versions:

40

mm

80

mm

- Standard Upholstery
- Premium Upholstery
- Omega Upholstery

Technical specification:

Upholstery: eco leather Highest quality and functionality

Liquid (water, blood) resistant High temp. resistant Soft, elastic and comfortable Scratch resistant Easy to clean

01 - Ecu

01 - Beige

02 - Blue

03 - Navy Blue

04 - Terracotta

05 - Brown

06 - Grey

07 - Graphite

08 - Black

09 - Red

11 - Honey

12 - Sea Green

13 - Orange

14 - Lime

15 - Light Blue

Technical specification:

Coating: 100% Polyurethane Phthalate - free

6 times more durable Low temp. resistant UV Resistant Chemical resistant Abrasion resistant

White - Dolce 6701

Ecu - Dolce 6703

Mango - Dolce 6731

Grass - Dolce 6722

Jet - Dolce 6724

Quetzal - Dolce 6734

Tangelo - Dolce 6718

Wine - Dolce 6738

Artesian - Dolce 6726

Graphite - Dolce 6739

Black - Dolce 6728

Chocolate - Dolce 6713

Technical specification:

Surface: 100% Polyurethane**Backing: 100% Polyester**

Resistance to hydrolis (8 years) No chemical solvents High scratch resistance Respects the environment
Durability and performance

OM006 - White

OM240 - Beige

OM315 - Sky

OM221 - Blue

OM314 - Cement

OM076 - Olive Grey

OM173 - Dark Brown

OM225 - Black

OM469 - Lima

OM157 - Forest

OM053 - Red

OM217 - Bordeaux

Accesories

Wedges

Bolsters

Paper holder

Protection mat

Stabilization Belt

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Type of table: Massage table

Number of sections: 2

Wheeling system: F4

Pivot: No

Product attributes:

Tables - Color: 00 - Ecu - Standard, 01 - Beige - Standard, 02 - Blue- Standard, 03 - Navy Blue - Standard, 04 - Teracota - Standard, 05 - Brown - Standard, 06 - Grey - Standard, 07 -

Length [cm]: 202
Width [cm]: 69
Height [cm]: 47-95
Weight [kg]: 95
Max. lifting capacity [kg]: 225
Max. working load [kg]: 300
Power supply [V / Hz]: 230 / 60
Warranty: 24 months
Medical device: Yes

Graphite - Standard, 08 - Black - Standard, 09 - Red - Standard, 11 - Honey - Standard, 12 - Marine - Standard, 13 - Orange - Standard, 14 - Lime - Standard, 15 - Navy - Standard, DOL-6701 White - Premium, DOL-6703 Ecru - Premium, DOL-6731 Mango - Premium, DOL-6722 Grass - Premium, DOL-6724 Jet - Premium, DOL-6734 Quetzal - Premium, DOL-6726 Artesian - Premium, DOL-6739 Graphite - Premium , DOL-6728 Black - Premium, DOL-6718 Tangelo - Premium, DOL-6738 Wine - Premium, DOL-6713 Chocolate - Premium
Tables - Upholstery thickness: 40 mm, 80 mm