

Prestige M-S2 - Massage table



Product short description:

Prestige M-S2 massage table is equipped with the innovative **Hallotronic®** system and Personal Authorization System for comfortable and convenient work with the patient.

The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with non-slip rubber feet that can be adjusted by up to 1 cm.

Prestige M-S2 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

Product gallery:

Product description:

Prestige M-S2 - Massage table

Prestige M-S2 massage and rehabilitation table is equipped with the innovative **Hallotronic®** system and **Personal Authorization System** for comfortable and convenient work with the patient. The headrest can be easily adjusted using a gas spring. It also features a contoured face opening with a plug. For levelling purposes, the table is equipped with non-slip rubber feet that can be adjusted by up to 1 cm.

Prestige M-S2 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

Prestige Series

A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Prestige Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

The Prestige M-S2 massage table features:

2-section massage and rehabilitation table

Adjustable head section from -70° to +40° with a gas spring.

Face hole in the head section with 2 removable face hole covers (round, crescent).

Adjustable leg section from -10° to 80° with 2 gas springs.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Personal Authorization Safety System with 2 safety magnetic keys.

HalloTronic® System for electric high/low height adjustment from 47 to 95 cm with a lifting capacity up to 225 kg

Types of base

You can choose between F0 and F4 wheeling system. F0 system composed of 2 legs and 2 small wheels allows for moving of Prestige Table by lifting it from one side. The system is recommended for users who relocate their Prestige Table less frequently. F4 system composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

F0 - Wheeling system

Composed of 2 legs with rubber feet and 2 small wheels for moving of Prestige Table by lifting it from one side. It is recommended for users who relocate their Prestige Table less frequently.

F4 - Wheeling system

Composed of 4 directional castors with central brake allows you for a very easy and fast relocation of Prestige Table with or without a patient.

Technologies

Personal Authorization Safety System

Tables have an integrated electronic control unit allowing for height adjustment and adjustment of the postural drainage option with an indication of stand-by mode (green light) and Personal Authorization Safety System (PASS). Each Prestige Table is equipped with two magnetic safety keys (PASS Keys) in order to prevent unauthorized use of a table which may lead to severe accidents and compensation claims.

Hallotronic® Technology

We care about your working comfort! That's why we invented the unique HalloTronic® System that enables high/low height adjustment via hallotrons without the use of mechanical switches. Thanks to HalloTronic® System you can choose the right height of your Prestige Table without effort and noise.

Prestige - Massage tables are widely used in:

Massage tables from the Prestige Series are a top-quality choice for sanatoriums, offering a comprehensive range of solutions designed to meet the needs of even the most demanding customers. Built with advanced technology and precision craftsmanship, these tables fully **support the recovery and relaxation process** of the client.

Upholstery in Prestige Series

Temperature resistant

Comfortable surface

Scratch resistant

Upholstery thickness available in two versions:

40

mm

80

mm

1. *Chlorophyll*
 2. *Organisms*
 3. *Organisms*
 4. *Organisms*
 5. *Organisms*
 6. *Organisms*
 7. *Organisms*
 8. *Organisms*
 9. *Organisms*
 10. *Organisms*
 11. *Organisms*
 12. *Organisms*
 13. *Organisms*
 14. *Organisms*
 15. *Organisms*
 16. *Organisms*
 17. *Organisms*
 18. *Organisms*
 19. *Organisms*
 20. *Organisms*
 21. *Organisms*
 22. *Organisms*
 23. *Organisms*
 24. *Organisms*
 25. *Organisms*
 26. *Organisms*
 27. *Organisms*
 28. *Organisms*
 29. *Organisms*
 30. *Organisms*
 31. *Organisms*
 32. *Organisms*
 33. *Organisms*
 34. *Organisms*
 35. *Organisms*
 36. *Organisms*
 37. *Organisms*
 38. *Organisms*
 39. *Organisms*
 40. *Organisms*
 41. *Organisms*
 42. *Organisms*
 43. *Organisms*
 44. *Organisms*
 45. *Organisms*
 46. *Organisms*
 47. *Organisms*
 48. *Organisms*
 49. *Organisms*
 50. *Organisms*
 51. *Organisms*
 52. *Organisms*
 53. *Organisms*
 54. *Organisms*
 55. *Organisms*
 56. *Organisms*
 57. *Organisms*
 58. *Organisms*
 59. *Organisms*
 60. *Organisms*
 61. *Organisms*
 62. *Organisms*
 63. *Organisms*
 64. *Organisms*
 65. *Organisms*
 66. *Organisms*
 67. *Organisms*
 68. *Organisms*
 69. *Organisms*
 70. *Organisms*
 71. *Organisms*
 72. *Organisms*
 73. *Organisms*
 74. *Organisms*
 75. *Organisms*
 76. *Organisms*
 77. *Organisms*
 78. *Organisms*
 79. *Organisms*
 80. *Organisms*
 81. *Organisms*
 82. *Organisms*
 83. *Organisms*
 84. *Organisms*
 85. *Organisms*
 86. *Organisms*
 87. *Organisms*
 88. *Organisms*
 89. *Organisms*
 90. *Organisms*
 91. *Organisms*
 92. *Organisms*
 93. *Organisms*
 94. *Organisms*
 95. *Organisms*
 96. *Organisms*
 97. *Organisms*
 98. *Organisms*
 99. *Organisms*
 100. *Organisms*

Accesorios

Wedges

Bolsters

Paper holder

Protection mat

Stabilization Belt

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Type of table: Massage table

Number of sections: 2

Wheeling system: F4

Pivot: No

Length [cm]: 202

Width [cm]: 69

Height [cm]: 47-95

Weight [kg]: 95

Max. lifting capacity [kg]: 225

Max. working load [kg]: 300

Power supply [V / Hz]: 230 / 60

Warranty: 24 months

Medical device: Yes

Product attributes: