

## Mixer type - B-130L, B-200L



### Product short description:

Mixer type B-130L, B-200L are used to prepare the mud mass for thermotherapy treatments. The device performs the following functions:

Automatic heating and maintaining the mud mass in the proper consistency at the set temperature

Automatic periodic mixing of the mud mass, which allows to maintain the proper consistency of the mud mass;

Pouring the mass from the tank through a special valve without the need to use a scoop.

### Product gallery:

### Product description:

## Mixer type - B-130L, B-200L

Mixer type **B-130L** and **B-200L** are used to prepare the mud mass for thermotherapy treatments. In addition, the control electronics automatically turns on the device after a preset time, which allows to get the prepared mass to start working by the operating staff and reduces the waiting time for the mass to be ready for treatment. The device is designed for use in sanatoriums, nursing homes and other care and health care facilities.

### Electronic controls allow programming of functions:

### Product codes:

Reference: -

EAN13: -

UPC: -

## Product features:

## Product attributes: