

Invacmed - Vacuum massage unit



Product short description:

This is a **vacuum massage unit** that comes with **4 vacuum electrodes**, each with a diameter of 60 mm. It has a pulsation frequency of 15-90 imp/min and a maximum vacuum of 400 mbar.

Product gallery:





Product description:

Invacmed - Vacuum massage unit

InVacMed - Vacuum Massage Unit is designed for vacuum massage and for attaching electrodes during electrotherapy treatments. The device allows combined electrotherapy and massage treatments in continuous mode with a constant vacuum value, as well as a treatment in pulsed mode with a variable vacuum adjustable up to 400 mbar. The device can be used for almost all conditions treatable by electrotherapy, facilitating the attachment of electrodes to the patient's body during procedures. In particular: reducing muscle tension in chronic neck pain and treating venous and lymphatic stasis.

Features:

2-channel universal vacuum therapy device that works with most of the apparatuses available on the market (including: Gymna, Chattanooga, EiE Otwock, Astar),

Screen set at an angle to facilitate working with the device,

Intuitive operation,

Possibility of applying current through suction cups with integrated electrodes,

Quick and practical attachment of self-sealing suction cups with electrodes in hard-to-reach places,

Independent connections for electrodes,

The ability to perform vacuum massage without application of currents to stimulate blood circulation and metabolism,

Smooth regulation of vacuum to adjust the mode of operation to meet the needs of the therapy,

Display showing the actual intensity level of the treatment,

Electronic work control and protection against excessive suction power,

Indication of water tank filling,

Power-saving operation

Compatibility:

With the purchase of the appropriate connection cable, **Invacmed** is compatible with devices from:

Astar,

Chattanooga,

EiE Otwock,

Gymna,

Standard accessories:

4 vacuum electrodes (60 mm) with viscose inserts,

Full set of electrode cables,

Power supply cable.

Optional accessories:

Electrodes (30 and 90 mm),

Viscose inserts (30 and 90 mm),

Electrotherapy device connecting cable.

Technical data:

Pulse frequency (pulsed mode) [pulses/min.]:	15 - 90
Vacuum [mbar]:	60 - 400
Number of electrodes:	4
Power supply [V/Hz]:	230 / 50-60
Current consumption [A]:	0,1
Dimensions (W x H x D) [cm]:	33,5 x 13,5 x 28,5
Weight: [kg]	3,4
Dehydrator reservoir capacity [ml]:	min. 80
Fuse [mA/T/V]:	125/250
Optical alarm level "H2O"	no more than 50 ml

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes: